

Distinguishing Pneumonia with X-ray of Chest via Deep Learning

made by Xiao-Ying Zhuang

Outline

Summary

Problem

Solution

Target

Flow Diagram



Training

Procedure

Data Source

Preprocessing

Training Models

Criteria

Comparison

Result

Application

Demo (flask)

Demo (tkinter)

Conclusion

Future Outlook

Reference



Problem

TIME

SPOTLIGHT STORY HOW QUANTUM COMPUTING WILL TRANSFORM OUR WORLD



The Coming Collapse of the U.S. Health Care System

IDEAS • HEALTH

The outbreak of the Covid-19 in late 2019 triggered a global epidemic of infectious diseases and paralyzed the global healthcare system.

When a patient with suspected Covid-19 is sent to the hospital, it is necessary to quickly judge its severity so that the amount of medical treatment can be properly allocated to prevent severe patients from running out of resources.

MAJOR ARTICLE • Rev. Soc. Bras. Med. Trop. 53 • 2020 • <https://doi.org/10.1590/0037-8682-0354-2020>

COPY

③ **Health system collapse 45 days after the detection of COVID-19 in Ceará, Northeast Brazil: a preliminary analysis**

Daniele Rocha Queiros Lemos [...] Luciano Pamplona de Góes Cavalcanti

ABOUT THE AUTHOR

The Economist

The global health-care collapse

January 16, 2023 • 9 min read

Solution




In practice, the diagnosis of pneumonia depends on the following things:


- ▲ Clinical symptoms (cough, sputum, fever, pain)
- ▲ Medical Images (Abnormal shadowing in the chest X-ray Image)

However, the interpretation of chest X-ray images will be more time-consuming. Besides, the more severe the lesion is, the easier it is to interpret. In other words, if the lesion is very small, there may be interpretation errors. And poor exposure imaging may also affect interpretation.


By using deep learning model (CNN), the problems mentioned above can be solved quickly and precisely. Because CNN model can “read” X-ray images pixel by pixel, which is almost impossible to human beings.




Target



Flow Diagram



Procedure



Label

- Normal
- Pneumonia
 - Virus
 - Bacteria

Preprocessing


- Rotation_range
- Width_shift_range
- Height_shift_range
- Zoom range
- resize

Training

- weight
- Epoch
- batch size
- EarlyStopping
- ModelCheckpoint



Data Source - from kaggle



Preprocessing

- rotation_rang = 30
- width_shift_range = 0.05~0.2
- height_shift_range = 0.05~0.2
- zoom_range = 0.2
- horizontal_flip = false

Training Models



- Weight = {0:3,1:1}
- Input_shape = 128~666
- Epoch = 30~45
- Batch size = 32~64
- Model_monitor = Recall
- ReduceLROnPlateau
 - Factor = 0.2~0.3
- Early_stopping =
 - monitor = loss
 - Min_delta = 1e-5~1e-2

Criteria




		Actual	
		True	False
Prediction	True	TP	FP
	False	FN	TN

- Accuracy = $(TP + TN) / (TP + FP + FN + TN)$
- Precision = $TP / (TP + FP)$
- Recall = $TP / (TP + FN)$
- F1 score = $2 / (\text{Precision}^{-1} + \text{Recall}^{-1})$



Training Models - Self-build CNN



```
model = Sequential()
model.add(Conv2D(32, (3, 3), input_shape = img.shape))
model.add(Activation("relu"))
model.add(MaxPool2D(pool_size=(2, 2)))

model.add(Conv2D(64, (6, 6)))
model.add(Activation("relu"))
model.add(MaxPool2D(pool_size=(2, 2)))


model.add(Conv2D(128, (9, 9)))
model.add(Activation("relu"))
model.add(MaxPool2D(pool_size=(2, 2)))

model.add(Flatten())
model.add(Dense(64))
model.add(Activation("relu"))
model.add(Dropout(0.5))
model.add(Dense(1))
model.add(Activation("sigmoid"))
```


Training Models - Keras applications

- ❖ Both ResNet50V2 and ResNet152V2 are versions of ResNet that have improved the vanishing gradient problem when dealing with deep networks.
- ❖ InceptionV3 combines multiple convolution operations (such as 1x1, 3x3, and 5x5 convolutions), having pooling layers of different sizes and types to improve model efficiency.
- ❖ MobileNet is a lightweight network architecture using depthwise separable convolutions that can run on mobile devices with limited hardware.
- ❖ EfficientNetV2S is a version of EfficientNet, which is mainly designed based on the combination of network depth, width and resolution.


Ensembling




Comparison - PR curve




inceptionV3




EfficientNetV2S



Ensembled




ResNet152V2




Self-build CNN

Comparison - best threshold finding




Result

Pneumonia
Normal




Result

Pneumonia
Normal





Result

Pneumonia
Normal



Data Source - from kaggle



Preprocessing

- rotation_range=15 ~ 30
- width_shift_range = 0.05~0.2
- height_shift_range = 0.05~0.2
- zoom_range = 0.1 ~ 0.4
- horizontal_flip = false

Training Models

- Weight = {0:1, 1:2}
- Input_shape = 224~299
- Epoch = 30~45
- Batch size = 32~64
- Model_monitor = accuracy
- EarlyStopping
 - monitor = loss
- ReduceLROnPlateau
 - factor = 0.2~0.3
 - min_delta = 1e-3~1e-5



Criteria




		Actual	
		True	False
Prediction	True	TP	FP
	False	FN	TN


- Accuracy = $(TP + TN) / (TP + FP + FN + TN)$
- Precision = $TP / (TP + FP)$
- Recall = $TP / (TP + FN)$
- F1 score = $2 / (\text{Precision}^{-1} + \text{Recall}^{-1})$




Comparison - PR curve




MobileNet



ResNet50



Ensembled





ResNet152V2

Self-build CNN


Result

Bacteria
Virus




Result

Bacteria
Virus



Result




Bacteria
Virus

Future Outlook



Process	Present Problem	The Way to Improve
Data Source	Resolution of X-ray imaging affects the accuracy of prediction.	Choose clear X-ray images.
Classification	Only divided into "Virus" & "Bacteria" when pneumonia is predicted.	Classify specific pathogens (e.g. tuberculosis) or virus type
Prediction	Unable to handle the coexistence of viruses and bacteria.	Use Normal / Bacteria & Normal / Virus to classify.
Ensembling	Extreme results affect the average.	Choose "vote" as the ensembling way but not "average".
Competitiveness	No other ready-made model to compare.	Compare with the result of the one generated by Teachable Machine.

Future Outlook



Reference

- ✓ <https://www.kaggle.com/datasets/paultimothymooney/chest-xray-pneumonia/code>
- ✓ <https://www.kaggle.com/code/karnikakapoor/pneumonia-diagnosis-convnet-model>
- ✓ <https://www.kaggle.com/code/madz2000/pneumonia-detection-using-cnn-92-6-accuracy>
- ✓ <https://www.kaggle.com/code/alisultanov/pneumonia-diagnosis-using-resne>
- ✓ <https://www.kaggle.com/competitions/vinbigdata-chest-xray-abnormalities-detection/data>
- ✓ <https://medlineplus.gov/xrays.html>
- ✓ <https://time.com/6246045/collapse-us-health-care-system/>
- ✓ <https://www.scielo.br/j/rsbmt/a/98LMbshKyrXVc7sC4nZwkSG/?lang=en>
- ✓ <https://finance.yahoo.com/news/global-health-care-collapse-190439531.html>